Precision Health Solutions

# Respiratory Panel

Nasopharyngeal

## SPECIMEN COLLECTION PROCEDURES

### SPECIMEN TYPE

- Nasopharyngeal swab submitted in Viral Transport Media (VTM) provided by Precision Health Solutions.
- Specimen should be submitted to the laboratory for testing as soon as possible, but no later than 2 days after collection.

### COLLECTION PROCEDURE

1. Wear personal protective equipment as required.
  - Gloves
  - Disposable lab coat
  - Eye protection
  - Droplet protection (Surgical mask or N95 mask if COVID-19 suspected)
2. Collect the NP swab.
  - Open the swab kit while ensuring sterility of swab.
  - Insert the swab into the nostril as far as you can. If you meet resistance, do not push through as you may cause damage to nasal mucosa.
  - Once inserted fully, hold in place for three seconds then slowly rotate for five seconds.
3. Withdraw the swab and place into the tube with media.
  - Ensure that you do not contaminate the swab while placing into the transport medium.
4. Break the swab at the break point.
5. Close the lid tightly, leaked specimens will not be processed.
6. Immediately label the sample collection container with first name, last name and date of birth.
  - *If this step is missed, your sample may be subject to rejection.*
7. Place the completed test requisition and any other supporting documentation in the pocket of the specimen bag.
8. Contact PHS to schedule courier pick-up.

Continue to the next page for detailed 'Oropharynx' collection instructions

info@precision-healthsolutions.com

727-235-0886


833-288-9397

Precision Health Solutions

## Respiratory Panel

Oropharynx

### SPECIMEN COLLECTION PROCEDURES

#### SPECIMEN TYPE

- Oropharynx swab submitted in Viral Transport Media (VTM) provided by Precision Health Solutions.
- Specimen should be submitted to the laboratory for testing as soon as possible, but no later than 2 days after collection.

#### COLLECTION PROCEDURE

1. Assemble all required supplies and equipment prior to collecting sample.
2. Position head slightly back, depress the tongue while visualizing the oropharynx.
3. Using the provided sterile swab, rub the swab up and down the back of the throat to collect cellular material where pathogens are at the highest concentration.
  - If multiple sites are present, collect sites and submit into one container.
  - Avoid touching the tongue and cheeks.
4. Place swab into VTM container.
5. Break the swab at the breakpoint.
6. Close the lid tightly, leaked specimens will not be processed.
7. Immediately label the sample collection container with first name, last name, and date of birth.
  - *If this step is missed, your sample may be subject to rejection.*
8. Place the completed test requisition and any other supporting documentation in the pocket of the specimen bag.
9. Contact PHS to schedule courier pick-up.

info@precision-healthsolutions.com

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